



Certificate of Conformance

ILS glass syringes are manufactured to exacting standards and are quality tested to meet or exceed our minimum performance specifications for accuracy and precision. Our testing methods are in accordance with GLP protocols and utilize equipment that is traceable to NIST standards. Copies of actual test results are available upon request. To maintain syringe accuracy always rinse with 5 volumes of water or methanol after use

Volume: 0,5µl, 1µl, 5µl, 10µl, 25µl, 50µl, 100µl, 250µl, 500µl, 1,0ml, 2,5ml, 5,0ml, 10ml

Accuracy: +/-1% of end volume

Stützerbach, 27.09.2011

Signature:.....

Heinz Ewald